

MANUAL

**TDCS FOR THE
TREATMENT OF
MAJOR DEPRES-
SIVE DISORDER:
FREQUENTLY
ASKED QUESTIONS**

SOOMA OY

Atomitie 5C

FI-00370 Helsinki, Finland

Tel. +358 10 328 9811

Email: info@soomamedical.com

www.soomamedical.com

MK40- 1.0 EN Sooma Oy 2022

CONTENTS

*	3	ABSTRACT
*	4	ABOUT SOOMA
*	8	USING TDCS IN THE TREATMENT OF MDD
*	10	THE TREATMENT
*	15	PATIENT SELECION CRITERIA
*	18	CONTRAINDICATIONS AND SIDE EFFECTS
*	20	USING SOOMA TDCS THERAPY
*	24	TESTIMONIALS
*	25	REFERENCES

QUESTIONS

1. What are the main differences between tDCS therapies and other neuromodulation treatment options for MDD?
2. How long does one course of treatment take?
3. How long does it take for the effects of the tDCS therapy to be noticed?
4. How long do the tDCS treatment effects last?
5. Can tDCS therapy be combined with other neuromodulation therapies?
6. How are other depressive symptoms such as anxiety or sleep disorders affected by tDCS therapy?
7. Who is a suitable patient for tDCS therapy?
8. How does tDCS therapy work in combination with other non-neuromodulation depression treatments?
9. Can patients with bipolar depression be treated with tDCS therapy?
10. Can tDCS be used to treat MDD in epileptic patients?
11. What are the contraindications for the use of tDCS in the treatment of MDD?
12. What are the side-effects associated with the use of tDCS therapy?
13. How long does it take to train the staff on using Sooma Depression Therapy?
14. How often does the staff need to be updated on how to use Sooma Depression Therapy?
15. How long will it take to train the patients on how to use Sooma Depression Therapy at home?
16. Is it safe for the patient to perform the treatment by themselves at home?
17. Is there a recommendation whether to start the treatment at the clinic or at home?
18. How can I ensure that my patients are following the prescribed treatment correctly?
19. How can I make sure that my patients return the device after the treatment period is over?

Keywords: *tDCS, MDD, depression, neuromodulation, treatment, therapy.*

ABOUT SOOMA

MISSION AND VISION

SOOMA OY IS A FINNISH MEDICAL DEVICE COMPANY ESTABLISHED IN 2013. SINCE THE START, WE HAVE BEEN WORKING CLOSELY WITH LEADING EXPERTS IN PSYCHIATRY AND CLINICAL NEUROPHYSIOLOGY TO DEVELOP EFFECTIVE TREATMENT SOLUTIONS FOR HEALTH CONDITIONS IN THESE AREAS.



TUOMAS NEUVONEN
CEO AND FOUNDER

Tuomas is a medical device entrepreneur with extensive expertise in brain stimulation devices. He is the mind behind the idea about the company and the product, and has built them up from the ground to a international business in over 30 countries.



MIKA NIKANDER
COO

Mika is our experienced expert with a portfolio of projects from product development, quality and regulatory affairs, and clinical support activities. He joined Sooma already in 2013 after finishing his own medical device product development project.



PETER HILGERT
CCO

With more than 20 years of experience in different commercial roles in the pharmaceutical industry, Peter has a profound knowledge of marketing and sales in large organizations, as well as extensive experience in the start-up environment.



SOOMA TDCS

In 2014, we launched Sooma Depression Therapy, a brain stimulation treatment for Major Depressive Disorder, along with Sooma tDCS, the high-end stimulation device that is used to conduct the treatment. Sooma tDCS medical devices are manufactured in Finland in accordance with the international ISO 13485 quality management system for medical devices. The treatment device developed by Sooma has received the following approvals from the world's health authorities, among others:

2014: Class IIA CE marking (Depression) | 2015: Malaysian Medical Device Authority (MDA) | 2016 Australian TGA (Therapeutic Goods Administration) | 2016: Health Canada (Chronic pain class II) | 2017: Singapore HSA (Health Science Authority) | 2017: CE marking (Chronic pain) | 2017: Turkey (Depression) | 2018: Mexican COFEPRIS (Federal Commission for Protection against Health Risks) | 2019: Indonesia (Depression)

TDCS FOR THE TREATMENT OF MAJOR DEPRESSIVE DISORDER: FREQUENTLY ASKED QUESTIONS

Transcranial direct current stimulation (tDCS) relieves the symptoms of major depressive disorder (MDD) by modulating cortical excitability with a weak current. tDCS delivers a constant current that induces changes in neuronal excitability in a polarity-dependent manner. In Sooma Depression Therapy, the positive anodal stimulation is used to increase neuronal excitability at the left dorsolateral prefrontal cortex (DLPFC), which is found to be hypoactive in depressed patients. The flow of current from the positive to the negative electrode influences the activity in the prefrontal cortex, relieving symptoms in depressed patients.

tDCS therapies are an increasingly used, safe, and effective method to treat MDD. However, in our daily practice we find the same doubts and questions regarding tDCS being raised time and again. The goal of this manual is to provide science-based answers to the most frequently asked questions regarding the use of tDCS for the treatment of MDD.



01

The treatment

In this section, we cover the most common questions regarding the use of tDCS therapies for the treatment of MDD.

We review the main differences between tDCS and other neuromodulation treatment options, as well as the possible synergic effects of those in combination with tDCS.

We cover the possible expected results of the treatment in terms of response to the treatment, duration of the response, and effect on common concomitant depressive symptoms such as anxiety and sleep disorders.

02

Patients selection criteria

In this section, we cover the most common questions regarding the criteria for determining which patients are the best suited to be treated with tDCS.

We describe the hypothetical ideal patient to be treated with tDCS and take a look at special cases, such as patients suffering from bipolar depression or other neurologic disorders.

We also present the feasibility of using tDCS in patients being treated with other non-neuromodulation treatments, and the possible synergic effects of such combination.

03

Contraindications and side-effects

In this section, we take a look at the contraindications and possible safety considerations regarding the use of tDCS therapies, with a special mention to the most common side-effects.

USING TDCS IN THE TREATMENT OF MDD

INTRODUCTION

TDCS RELIEVES THE SYMPTOMS OF MDD BY MODULATING CORTICAL EXCITABILITY OF THE DORSOLATERAL PREFRONTAL CORTEX WITH A WEAK STIMULATION CURRENT.

The use of neuromodulation techniques for the treatment of mood disorders has been ongoing since the beginning of the 20th century.

Starting with Electroconvulsive Therapy (ECT) and moving towards a wider range of invasive and non-invasive therapies, the use of electricity or electromagnetic fields to regulate neuronal activity is a proven method to treat diverse psychiatric conditions.

Invasive neuromodulation involves surgically implanting several components into the body: a small stimulation device –often placed in the chest area; electrodes, placed in the stimulation targets; and leads connecting the two. The stimulation device sends stimulus, in the form of electric impulses, through the leads to the electrodes, stimulating the target area for the desired result. Examples of invasive neuromodulation are Deep Brain Stimulation (DBS), Spinal Cord Stimulation (SCS), and Vagus Nerve Stimulation (VNS), which are most commonly used to treat movement disorders, pain, and epilepsy, respectively.

In **non-invasive neuromodulation**, the stimulation electrodes or the electromagnetic coil are placed externally on top of the scalp, in such locations that the electric field will be formed to the target areas of the brain. To treat major depressive disorder (MDD), ECT and Transcranial magnetic stimulation (TMS) are currently the most commonly used methods. Transcranial direct current stimulation (tDCS) is a relatively newer stimulation method to treat depression and is the method used by Sooma Therapies.

tDCS relieves the symptoms of MDD by modulating cortical excitability with a weak current. tDCS delivers a constant current that induces changes in neuronal excitability in a polarity-dependent manner (for review see Tortella et al. 2015). In Sooma Depression Therapy, the positive anodal stimulation is used to increase neuronal excitability at the left dorsolateral prefrontal cortex (DLPFC), which is found to be hypoactive



tDCS is an effective and well-tolerated treatment for MDD with an A level of evidence

	Condition	Recommendation	Boiled effect sizes*
Pain	Neurogenic pain	Avoid tDCS probably effective in reducing neurogenic pain (Level B)	-0.29 [-0.40, 0.02]
	Migraine	Avoid tDCS probably effective in reducing migraine pain (Level B)	-0.62 [-1.23, -0.01]
	Migraine	Avoid tDCS probably effective in reducing migraine pain (Level B)	-0.41 [-1.40, 0.59]
MFS	Postoperative acute pain	No recommendation	Not estimable
	Postoperative acute pain	Postoperative avoid tDCS probably effective in reducing postoperative analgesia and pain (Level B)	-0.70 [-1.00, -0.30]
	Low back pain	No recommendation	Not estimable
FD	Motor function	Avoid tDCS of ipsilateral M1 possibly effective for motor function in PD (Level C); avoid prefrontal tDCS probably not effective for motor function in TD (Level B)	-0.88 [-0.88, -0.88]
	Cognitive function	Avoid tDCS probably effective for cognitive function in PD (Level B)	-0.33 [-1.00, 0.35]
	Chronic stroke	Motor function	Avoid tDCS of ipsilateral M1 probably effective for motor rehab in chronic stroke (Level B)
Aphasia	Motor function	Avoid tDCS of ipsilateral M1 to enhance robotic therapy probably not effective for motor rehab in chronic stroke (Level B)	0.23 [-0.31, 0.77]
	Aphasia	Cathodal tDCS of contralateral M1 probably effective for motor rehab in chronic stroke (Level B)	0.41 [-0.18, 1.00]
	Aphasia	Bilateral tDCS of M1 probably effective for motor rehab in chronic stroke (Level B)	0.44 [0.15, 0.75]
Subacute stroke	Motor function	Avoid tDCS of Broca's area, its homologues, or Wernicke's area is possibly effective in chronic post-stroke aphasia rehab (Level C)	Broca's Area 0.43 [0.23, 1.03]
	Motor function	Cathodal right frontotemporal tDCS is possibly effective in chronic post-stroke aphasia rehab (Level C)	Not estimable
	Motor function	Bilateral tDCS of Broca's area and cathodal stimulation of its homologues is possibly effective in chronic post-stroke aphasia rehab (Level C)	Not estimable
Epilepsy	Motor function	Avoid tDCS of ipsilateral M1 probably effective for motor rehab in subacute stroke (Level B)	0.44 [-0.03, 0.91]
	Motor function	Avoid tDCS of ipsilateral M1 to enhance robotic therapy is definitely not effective for motor rehab in subacute stroke (Level A)	0.51 [-0.33, 0.31]
	Motor function	Cathodal tDCS of contralateral M1 probably effective for motor rehab in subacute stroke (Level B)	0.47 [0.18, 0.86]
MDD	Depression	Avoid tDCS probably effective for motor rehab in subacute stroke (Level B)	Not estimable
	Depression	Avoid tDCS of M1 is possibly effective for motor rehab in subacute stroke (Level C)	0.89 [-0.07, 0.81]
	Depression	Bilateral tDCS of M1 is possibly effective for motor rehab in subacute stroke (Level C)	0.89 [-0.07, 0.81]
OCD and TS	OCD	Cathodal tDCS probably safe (no increase in seizures) and effective (decrease in seizures) in epilepsy (Level B)	-0.70 [-1.38, -0.02]
	OCD	Avoid tDCS probably effective in improving OCD symptoms (Level C)	-0.41 [-0.66, -0.06]
Schizophrenia	TS	No recommendation	Not estimable
	Schizophrenia	Avoid tDCS of Broca's area, its homologues, or Wernicke's area is possibly effective in chronic post-stroke aphasia rehab (Level C)	-0.47 [-0.72, -0.22]
Addiction	Alcohol craving or relapse	Right tDCS avoid or left tDCS avoid probably effective in decreasing relapse or craving in alcohol addiction (Level B)	0.21 [-0.68, 0.49]
	Crack-cocaine	No recommendation	Not estimable
Addiction	Crack-cocaine	No recommendation	Not estimable
	Methylphenidate	No recommendation	Not estimable

Summary of Recommendations on tDCS Efficacy According to Clinical Indications (Fregni et al. 2021)

in depressed patients (Fitzgerald et al. 2006, Grimm et al. 2008). The flow of current from the positive to the negative electrode influences the neuronal activity in the prefrontal cortex (Tortella et al. 2015), and leads to relieving symptoms in depressed patients (Brunoni et al. 2010).

tDCS can be administered as a monotherapy, or as an adjunct treatment to enhance the effect of pharmaceutical or psychological therapy (Brunoni et al. 2013, Segrave et al. 2014). The efficacy and safety of tDCS have been shown in randomized controlled clinical studies, either as a monotherapy or as an adjunct therapy to pharmaceutical or psychosocial treatments (Brunoni et al. 2013, National Institute for Health and Care Excellence 2015a, National Institute for Health and Care Excellence 2015b, Brunoni et al. 2011, Mutz et al. 2018, Brunoni et al. 2016).

Based on evidence from clinical trials and meta-analyses, tDCS is an effective and well-tolerated treatment for MDD with an A level of evidence (Fregni et al. 2021).

THE TREATMENT


BASICS OF TDCS

WHAT ARE THE MAIN DIFFERENCES BETWEEN TDCS THERAPIES AND OTHER NEUROMODULATION TREATMENT OPTIONS FOR MDD?

tDCS therapies are a non-invasive neuromodulation method that uses a weak electric current (approx. 2mA) to stimulate targeted areas of the brain. The treatment protocol often involves a daily stimulation session for 3-6 weeks, which totals about 15-30 sessions. tDCS is regarded as a safe treatment method and does not cause any serious adverse effects (Brunoni et al. 2010, Segrave et al. 2014). Non-serious side effects include mild headache, and an itching sensation or redness in the stimulated area (Russo et al. 2017). In addition to using tDCS as monotherapy or augmentation therapy, tDCS can be used as a maintenance treatment to retain the effects after acute treatment with either Electroconvulsive Therapy (ECT) or Transcranial Magnetic Stimulation (TMS). More importantly, certain tDCS therapies—like Sooma Therapy, can be administered at home by the patient, alleviating the workload of clinics and hospitals.

In contrast, other non-invasive therapies like ECT or TMS, require the patient to visit the hospital or clinic daily for the administration of the treatment either because the patient needs to be under medical supervision due to anesthesia (ECT), risk of seizures (TMS), and just simply because the equipment is not mobile and the treatment is not suitable for self-administration. Also, associated risks and side-effects, even when uncommon, are more serious (Taylor, Galvez, and Loo 2018, Andrade, Arumugham & Thirthalli 2016).

Finally, invasive treatments like Deep brain stimulation (DBS), Spinal Cord Stimulation (SCS), and Vagus Nerve Stimulation (VNS) require a surgical procedure during which implants are inserted into the patient's body, and they carry the risks associated with demanding invasive surgery and post-surgery recovery as well as a wider range of associated side-effects. Furthermore, the first two treatments require subsequent surgeries for the replacement of the device batteries every 2-10 years, depending on the device and battery type.




VNS implant on a patient's chest (Image: Shutterstock)

Neuromodulation technique	Treatment protocol	Location of session	Anesthesia required	Hospitalization required	Side-effects and associated risks	Follow ups
Invasive neuromodulation						
Deep Brain Stimulation (DBS)	Pre-surgery evaluations and test. Surgical implant installation with a few days hospitalization. Programming of the device a few weeks after surgery. Adjustment of the protocol over a few weeks to months.	Hospital (several visits before, during and after implant surgery) and home (stimulation)	Yes	Yes	Risks associated with major surgery and general anesthesia, seizure, headache, confusion, difficulty concentrating, hardware complications.	Need to change the device batteries every 3-5 years (new surgery required)
Vagus Nerve Stimulation (VNS)	Pre-surgery evaluations and test. Surgical implant installation with a few days hospitalization. Programming of the device a few weeks after surgery. Adjustment of the protocol over a few weeks to months.	Hospital (several visits before, during and after implant surgery) and home (stimulation)	Yes	Yes	Risks associated with major surgery and general anesthesia, difficulty swallowing, vocal cord paralysis, throat pain/hoarseness, tingling or prickling of the skin, shortness of breath.	Need to change the device batteries every 3-8 years (new surgery required)
Non-invasive neuromodulation						
Electro-Convulsive Therapy (ECT)	Pre-treatment evaluation. 1 hour treatment sessions 2/3 times a week during 3/4 weeks.	Hospital	Yes	Temporary (outpatient)	Risks associated with seizure induction and anesthesia, loss of memory, headaches, nausea, muscle pain or soreness, disorientation and confusion.	Maintenance sessions either by ECT, TMS or tDCS therapies.
Transcranial Magnetic Stimulation (TMS)	Pre-treatment evaluation. 1 hour treatment, 5 times a week, during 4/6 weeks.	Hospital	No	Temporary (outpatient)	Risk of seizures, maniac switch, hearing loss. Headaches, skin discomfort during stimulation, muscle spasms, lightheadedness.	Maintenance sessions either by TMS or tDCS therapies.
Transcranial Direct Current Stimulation (tDCS)	30 minutes treatment, 5 times a week, for a minimum of 3 weeks.	Hospital or home	No	No	Redness of the stimulation area, itching or tingling sensation on the skin, transient headache.	Maintenance sessions by tDCS if needed.

HOW LONG DOES ONE COURSE OF TREATMENT TAKE?

One course of treatment usually lasts a minimum of 3 weeks, depending on the patient profile, after which maintenance treatment can be applied. Generally, a longer treatment period—that is a higher number of sessions, correlates with a better response.



According to our clinical data, the minimum recommended treatment course with Sooma Depression Therapy is 3 weeks of acute treatment, followed up by a maintenance period. However, the clinical practice of most of our customers seem to indicate that a longer period treatment will yield better results. During the acute treatment phase, a 30-minutes stimulation session is performed five times a week.

HOW LONG DOES IT TAKE FOR THE EFFECTS OF THE TDCS THERAPY TO BE NOTICED?

The successful response to the therapy will depend on the characteristics of the patient such as the severity of the depressive episode, the type of depressive episode (unipolar versus bipolar), or a concurring treatment with antidepressants or other drugs (Brunoni et al. 2016). Although some patients can receive significant improvement in their condition already starting after one week of tDCS therapy, in most cases the response will require several weeks of treatment (Gorrigk et al., 2020). Furthermore, the level of response is correlated with tDCS dose (the number of sessions, their duration and current level used), and it is inversely associated with existing treatment-resistance (Brunoni et al. 2016).

Generally speaking, the treatment-resistant MDD patients are a challenging patient group also for tDCS therapy, and the three clinical studies using a one, two or three weeks tDCS treatment protocol, respectively, each failed to find significant improvement with tDCS therapy in the treatment-resistant MDD patients (Bennabi et al, 2014; Palm et al., 2011 and Blumberger et al., 2012). It was concluded after these studies that more intensive (more sessions, longer treatment period, etc) tDCS treatment could be required in this chronic and harder to treat depression group.

However, it must be considered that the response rate is very patient-specific and that the more adequate the patient selection process, the faster and more optimal the results of the therapy will be. More on the patient's selection criteria will be explored in upcoming sections.



HOW LONG DO THE TDCS TREATMENT EFFECTS LAST?

Regarding the duration of therapeutic effects on the long-term, there is scarce literature with follow ups after acute treatment sessions for MDD patients. Valiengo et al. (2013) analyzed the 24-weeks follow-up period after a two-week-acute tDCS treatment (10 sessions) followed by maintenance treatment with one tDCS session every other week. They found that 23% of the nontreatment-resistant and 90% of the treatment-resistant MDD patients relapsed during the follow-up, the mean response duration being 11,7 weeks (Valiengo et al, 2013).

Later, Aparicio et al. (2019) did a 6-month's follow up study where the patients received three weeks of acute tDCS treatment (15 sessions), followed by maintenance treatment with two tDCS sessions every week for relapse prevention. Similarly to previous study, they found a lower relapse trend for the nontreatment-resistant vs. the treatment-resistant patients, the rates of relapse being 7.7% and 45.5%, respectively. The mean response duration in their study was 17,5 weeks.

Interestingly, these data suggest that longer acute tDCS schedules with more intense maintenance treatment improve prevention of relapse and underline the importance of scheduling long enough tDCS treatment periods especially to the more severe and treatment-resistant patients.



CAN TDCS THERAPY BE COMBINED WITH OTHER NEUROMODULATION THERAPIES?

Absolutely. There is evidence both in the literature and clinical practice that patients that respond positively to other non-invasive neuromodulation treatments will also respond well to tDCS as a maintenance treatment option (Cha et al. 2016). Furthermore, tDCS can be used as preconditioning treatment prior to the administration of an acute phase treatment with rTMS, as it has been demonstrated that tDCS preconditioning enhances cortical plasticity and can shape the direction of rTMS-induced after-effects. (Lang et al. 2004, Consentino et al. 2012, Alkhasli et al. 2022).

HOW ARE OTHER DEPRESSIVE SYMPTOMS SUCH AS ANXIETY OR SLEEP DISORDERS AFFECTED BY TDCS THERAPY?

In a recent publication, Goerigk et al (2021) analysed the therapeutic effect to four different depression symptoms clusters: core depressive symptoms, atypical, sleep/insomnia and guilt/anxiety. Among the symptoms' clusters, tDCS treatment results were superior to placebo in core depressive and sleep/insomnia scores. On the other hand, tDCS and pharmacotherapy seemed to differ to some extent in their efficacies against different symptoms of depression. tDCS was effectively reducing the core depressive and sleep/insomnia cluster symptoms whereas escitalopram was reducing the core depressive and guilt/anxiety cluster symptoms. Also the efficacy of tDCS on the guilt/anxiety cluster symptoms was better than with placebo, but due to high variability this difference did not reach statistical significance.

Sleep/ insomnia	Atypical	Core depressive	Guilt/ Anxiety
Sleep-onset insomnia	Suicide	Mood (sad)	Guilt and delusions
Early morning insomnia	Psychomotor agitation	Loss of interest	Energy/ Fatigability
Midnocturnal insomnia	Psychomotor slowing		Reduced libido
	Hypochondriasis		Psychological anxiety
			Somatic anxiety

MDD symptoms clusters as per HAM-scale.

PATIENT SELECTION CRITERIA

USUAL AND SPECIAL CASES

ANY PATIENT SUFFERING FROM MDD WHO IS NOT OTHERWISE INCLUDED IN THE CONTRAINDICATIONS SECTION OF THIS DOCUMENT CAN BE TREATED WITH TDCS THERAPY

WHO IS A SUITABLE PATIENT FOR TDCS THERAPY?

As a general rule, any patient suffering from MDD who is not otherwise included in the contraindications section of this document can be treated with tDCS therapy. However, it is worth noting that the patient's own specifications will have an effect on the final result of the treatment.

According to the literature, tDCS has lesser efficacy in treatment-resistant depressed patients (Bennabi & Haffen 2018, Brunoni et al. 2016). On the other hand, depression severity is positively related to clinical improvement (Bennabi & Haffen 2018).

Therefore, it could be stated that tDCS is suitable for patients with a mild-to-severe form of MDD who do not meet the criteria for resistant depression (Bennabi & Haffen 2018).

HOW DOES TDCS THERAPY WORK IN COMBINATION WITH OTHER NON-NEUROMODULATION DEPRESSION TREATMENTS?

The combination of treatment methods has been researched and showed that the final effects of tDCS are affected to a great extent by using concomitant pharmacotherapy and/or other interventions, such as cognitive control therapy (CCT).

Regarding the use of antidepressant medication, setraline has been shown to have a synergistic therapeutic effect when combined with tDCS (Brunoni et al. 2013a). On the other hand, the use of benzodiazepines has been correlated with a worse improvement through neuromodulation treatments such as tDCS and rTMS, which could be due to the fact that cortical excitability is decreased when using benzodiazepines (Goerigk et al. 2020).

Concerning the effects of combining tDCS with cognitive therapy, the results are mixed. While Segrave et al. (2014) found preliminary evidence that concurrent CCT enhances antidepressant outcomes from tDCS, Brunoni et al. (2014) could not find significant relevance on the superiority of combined active treatments versus sham.

CAN TDCS BE USED TO TREAT MDD IN EPILEPTIC PATIENTS?

The safety of tDCS has been extensively studied in regards to epilepsy across multitude of stimulation paradigms. There are no reports of epileptic seizures caused by tDCS in the literature (Bikson et al. 2016). tDCS has also been used for epileptic patients. In a meta-analysis of 6 articles including 65 stimulated subjects, San Juan et al. (2015) found evidence that tDCS are shown to be preliminary safe and effective for patients with epilepsy.

Furthermore, the evidence points towards cathodal tDCS leading to reduction in seizure frequency (Regner et al. 2018) and may lead to better seizure control in drug-resistant focal epilepsy (Sudbrack-Oliveira et al. 2021). Therefore, Sudbrack-Oliveira et al. (2021) concluded that cathodal tDCS is overall safe and not directly associated with seizures in adults and children with drug-resistant epilepsy.



CAN PATIENTS WITH BIPOLAR DEPRESSION BE TREATED WITH TDCS THERAPY?

Bipolar depression is challenging to maintain in balance and unfortunately the available antidepressant treatment options are unsatisfactory and carry the risk of switching from depression to hypomania/mania (Roger et al. 2019). Thus, there is an apparent need to develop new antidepressant treatments for patients with bipolar depression.

tDCS treatments have shown promise as an effective and safe treatment option for bipolar depression (Sampaio-Junior et al. 2018, Dondé et al. 2017, Dondé et al. 2018). However, as with other antidepressant treatments, a risk for hypomania/mania switch has been reported with tDCS treatment of bipolar depression, suggesting that concomitant use of mood stabilizers and closer monitoring of the patients is prudent (Dondé et al 2017, Dondé et al. 2018).

Regarding the effects of tDCS stimulation in HDRS-17 Score on bipolar patients, Sampaio-Junior et al. (2018) demonstrated that, over the course of 6 weeks, the active tDCS stimulation was superior to sham. Active stimulation provided a steady decrease from the baseline score, that in week 6 was more than 10 points lower than in the beginning of the research. On the other hand, while sham stimulation initially also caused a decrease of the score, this was smaller, and was not sustained over the duration of the experiment, increasing again after the fourth week. Furthermore, active tDCS was superior to sham for providing sustained response (reduction in depression score of more than 50%), showing that the effect of active tDCS was also clinically significant in this patient group.

Most importantly, some emerging data suggest that tDCS improves neurological soft signs and sleep quality in euthymia, and therefore can be useful for relapse prevention (Dondé et al. 2018).

CONTRAINDICATIONS

AND SIDE-EFFECTS

WHAT ARE THE CONTRAINDICATIONS FOR THE USE OF TDCS IN THE TREATMENT OF MDD?

Most patients are eligible for tDCS therapy, but the eligibility of those suffering from some very specific comorbidities and acute conditions must be consulted with a specialist.

The list of widely accepted contraindications includes patients who present:

- pacemakers

Although scientific studies have shown that a mild electrical current administered on the scalp will not affect a pacemaker (Roncero et al. 2020), it is a vital implant to the patient, and thus extra care is necessary.

- metal inside the skull

An object inside or at the skull is generally a contraindication for receiving electrical currents, since there is a risk that those will damage the brain tissue surrounding the metallic implant. If the metal is located in the neck or the jaw, for instance dental fillings, it is safe to use tDCS.

- Acute eczema on the stimulation area


Broken skin on the stimulation area is a contraindication for using the treatment, as the current, even if mild, can worsen the condition.

WHAT ARE THE SIDE-EFFECTS ASSOCIATED WITH THE USE OF TDCS THERAPY?

No serious side effects or long term effects have been associated with the use of tDCS treatments. However, minor side effects during the stimulation might include (Liu et al. 2017):

1. Burning or tingling sensation on the electrode placement area
2. Transient headache
3. Redness of the scalp

The real-life effectiveness and tolerability of Sooma tDCS™ Depression Therapy was evaluated in a cohort of 302 MDD patients collected from ten primary and secondary care clinics worldwide. All patients were treated according to the standard Sooma tDCS™ stimulation protocol. **There were no serious adverse events** in the cohort during the 4,022 treatment sessions. The majority of reported adverse events were mild and transient reactions to the treatment.



USING SOOMA TDCS THERAPY FOR TREATING MDD

PRACTICALITIES

A TYPICAL TRAINING PROGRAM FOR NEW USERS OF SOOMA DEPRESSION THERAPY CONSISTS OF TWO ONE-HOUR SESSIONS. TYPICALLY, THERE IS NO NEED TO UPDATE THE TRAINING AFTERWARDS.

HOW LONG DOES IT TAKE TO TRAIN THE STAFF ON USING SOOMA DEPRESSION THERAPY?

A typical training program for people new to Sooma Depression Therapy consists of two one-hour sessions. The first session provides a general overview on what tDCS is, how it can help patients, along with a very brief overview of the scientific background Sooma treatment is based on. The second session consists of a workshop, where nurses and doctors get hands-on practice in administering the treatment. We arrange training sessions for medical staff on-demand. You can [contact us](#) to arrange a training session at your convenience.

HOW OFTEN DOES THE STAFF NEED TO BE UPDATED ON HOW TO USE SOOMA DEPRESSION THERAPY?

Generally, it is not necessary to update the staff about the device. However, if some new uses or changes occur, we will publish new training videos on our [online learning platform](#) for you to access them easily and at your convenience.

HOW LONG WILL IT TAKE TO TRAIN THE PATIENT ON HOW TO USE SOOMA DEPRESSION THERAPY AT HOME?

Typically, a one-hour-long session is enough for the patient to get a correct understanding on how the treatment is administered. However, up to 3 sessions can be



done onsite in the clinic for the patient to be confident enough with the device, after which the treatment can be continued at home.

If the patients have any doubts during the self-administration of the treatment, they can access our [Patient Material Bank](#) at any time to check how-to videos and printed materials.

IS IT SAFE FOR THE PATIENT TO PERFORM THE TREATMENT BY THEMSELVES AT HOME?

Absolutely. Given the absence of serious side-effects, tDCS therapies can be safely self-administered by anyone who has previously received adequate teaching by a medical professional, and provided that the treatment has been prescribed by a doctor or specialist who will oversee the establishment of an ongoing communication, frequent feedback, and controlled device settings (Alonzo et al. 2019).


IS THERE A RECOMMENDATION WHETHER TO START THE TREATMENT AT THE CLINIC OR AT HOME?

Although the goal of Sooma DEpression Therapy is to decongest the clinics by providing a treatment solution that can be easily used at home by the patients, it is recommended to have at least one first session under the supervision of a health professional. Doing the first session at the clinic can help to ensure that the patient (or a caregiver) knows how to use the device and is a good candidate for the treatment.

Typically, a one-hour-long session is enough for the patient to get a correct understanding on how the treatment is administered. However, it can be done for up to 3 sessions onsite in the clinic for the patient to be confident enough with the device, after which the treatment can be continued at home.

HOW CAN I ENSURE THAT MY PATIENTS ARE FOLLOWING THE PRESCRIBED TREATMENT CORRECTLY?

The ComfoPads that we provide are a reliable means of monitoring compliance: as they are single use and you provide the patient with the exact amount of disposables that they would need for a course of treatment, you can make sure that they will only perform the amount of sessions prescribed.



HOW CAN I MAKE SURE THAT MY PATIENTS RETURN THE DEVICE AFTER THE TREATMENT PERIOD IS OVER?

We have a loan agreement document in place that you can give for the patient to sign prior to starting the treatment. Furthermore, we recommend charging the patient a small deposit fee for the rental of the device which is adequate to the price level of the market, but still taxing to the patient. Such deposit will then be reimbursed upon the end of the treatment if the device is returned in a good condition.



TESTIMONIALS

THESE TESTIMONIALS HAVE BEEN GIVEN BY MEDICAL PROFESSIONALS ALL AROUND THE WORLD WHO USE SOOMA TDCS TO TREAT PATIENTS WITH MDD



Anu Kinnunen, MD and Chief Psychiatrist

North Karelia Hospital District (Siun Sote), Finland

“In our experience, tDCS has also enhanced the recovery of the hospitalized depressed patients - and if the patient has benefited from tDCS treatment, the treatment initiated in the hospital has been continued at home after discharge”

Margus Lõokene, MD and Senior Psychiatrist

North Estonia Medical Centre, Estonia

“The first thing we notice is the reduced anxiety [...] later on patients get more emotional, active, and after treatment we often get feedback from patients that they are functioning better in every day life, including work and social activities”

Rechdi Ahdab MD, PhD, Associated Professor of Medicine

Lebanese American University, Lebanon

“I have had an excellent experience with Sooma tDCS therapy specifically. The device is user friendly, easy to operate and reliable. Given its excellent safety profile, I mention this therapeutic options to my patients”

REFERENCES

LITERATURE CITED

- Alkhasli, I., Mottaghy, F. M., Binkofski F., et al. (2022) Preconditioning Prefrontal Connectivity Using Transcranial Direct Current Stimulation and Transcranial Magnetic Stimulation, Research Square.
- Alonzo, A., Fong, J., Ball, N., Martin, D., Chand, N., & Loo, C. (2019). Pilot trial of home-administered transcranial direct current stimulation for the treatment of depression. *Journal of affective disorders*, 252, 475–483.
- Andrade, C., Arumugham, S. S., & Thirthalli, J. (2016). Adverse Effects of Electroconvulsive Therapy. *The Psychiatric clinics of North America*, 39, 3, 513–530.
- Aparicio, L. V. M., Rosa, V., Razza, L. M., et al. Transcranial direct current stimulation (tDCS) for preventing major depressive disorder relapse: Results of a 6-month follow-up. *Depress Anxiety*. 2019, 36, 262– 268.
- Bennabi, D., Nicolier, M., Monnin, J., Tio, G., Pazart, L., Vandell, P., Haffen, E. (2014) Pilot study of feasibility of the effect of treatment with tDCS in patients suffering from treatment-resistant depression treated with escitalopram. *Clin Neurophysiol*, 126, 6, 1185-1189.
- Bennabi, D., & Haffen, E. (2018). Transcranial Direct Current Stimulation (tDCS): A Promising Treatment for Major Depressive Disorder?. *Brain sciences*, 8, 5, 81.
- Bikson, M., Grossman, P., Thomas, C., Zannou, A. L., Jiang, J., Adnan, T., Mourdoukoutas, A. P., Kronberg, G., Truong, D., Boggio, P., Brunoni, A. R., Charvet, L., Fregni, F., Fritsch, B., Gillick, B., Hamilton, R. H., Hampstead, B. M., Jankord, R., Kirton, A., Knotkova, H. ,... Woods, A. J. (2016). Safety of Transcranial Direct Current Stimulation: Evidence Based Update 2016. *Brain stimulation*, 9, 5, 641–661.
- Blumberger, D. M., Tran, L. C., Fitzgerald, P. B., Hoy, K. E., Daskalakis, Z. J. (2012) A randomized double-blind sham-controlled study of transcranial direct current stimulation for treatment-resistant major depression. *Front Psychiatry*. 17, 3, 74.
- Brunoni, A. et al. (2010) Neuromodulation approaches for the treatment of major depression: challenges and recommendations from a working group meeting. *Arq. Neuropsiquiatr*. 68, 433–451.
- Brunoni, A. et al. (2011) A systematic review on reporting and assessment of adverse effects associated with transcranial direct current stimulation. *Int. J. Neuropsychopharmacol*. 14, 1133–1145.
- Brunoni, A. et al. (2013) The Sertraline vs Electrical Current Therapy for Treating Depression Clinical Study. *JAMA Psychiatry* 70, 383.
- Brunoni, A. R. et al. (2016) Transcranial direct current stimulation for acute major depressive episodes: meta-analysis of individual patient data. *The British Journal of Psychiatry*, 208, 6, 522–531.
- Brunoni, A. R., Boggio, P. S., De Raedt, R., Benseñor, I. M., Lotufo, P. A., Namur, V., Valiengo, L. C., & Vanderhasselt, M. A. (2014). Cognitive control therapy and transcranial direct current stimulation for depression: a randomized, double-blinded, controlled trial. *Journal of affective disorders*, 162, 43–49.

- Cha, Y., Urbano, D., & Pariseau, N. (2016) Randomized Single Blind Sham Controlled Trial of Adjunctive Home-Based tDCS after rTMS for Mal De Debarquement Syndrome: Safety, Efficacy, and Participant Satisfaction Assessment. *Brain Stimulation*, 9, 4, 537-544.
- Cosentino, G., Fierro, B., Paladino, P., Talamanca, S., Vigneri, S., Palermo, A., Giglia, G., & Brighina, F. (2012) Transcranial direct current stimulation preconditioning modulates the effect of high-frequency repetitive transcranial magnetic stimulation in the human motor cortex. *European Journal of Neuroscience*, 35, 119–124.
- Dondé, C., Amad, A., Nieto, I., Brunoni, A. R., Neufeld, N. H., Bellivier, F., Poulet, E., & Geoffroy, P. A. (2017) Transcranial direct-current stimulation (tDCS) for bipolar depression: A systematic review and meta-analysis, *Progress in Neuro-Psychopharmacology and Biological Psychiatry*, 78, 123-131.
- Dondé, C., Neufeld, N.H. & Geoffroy, P.A. (2018) The Impact of Transcranial Direct Current Stimulation (tDCS) on Bipolar Depression, Mania, and Euthymia: a Systematic Review of Preliminary Data. *Psychiatr*, 89, 855–867.
- Fitzgerald, P. et al. (2006) An analysis of functional neuroimaging studies of dorsolateral prefrontal cortical activity in depression. *Psychiatry Res*. 148, 33–45.
- Fregni, F. et al. (2021) Evidence-Based Guidelines and Secondary Meta-Analysis for the Use of Transcranial Direct Current Stimulation in Neurological and Psychiatric Disorders, *International Journal of Neuropsychopharmacology*, 24, 4, 256–313.
- Goerigk, S. A., Padberg, F., Bühner, M., Sarubin, N., Kaster, T. S., Daskalakis, Z. J., Blumberger, D. M., Borriore, L., Razza, L. B., Brunoni, A. R. (2020) Distinct trajectories of response to prefrontal tDCS in major depression: results from a 3-arm randomized controlled trial. *Neuropsychopharmacology*. 46, 4, 774-782.
- Goerigk, S. A., Padberg, F., Chekroud, A., Kambeitz, J., Bühner, M., & Brunoni, A. R. (2021). Parsing the antidepressant effects of non-invasive brain stimulation and pharmacotherapy: A symptom clustering approach on ELECT-TDCS, *Brain Stimulation*, 14, 4, 906-912.
- Grimm, S. et al. (2008) Imbalance between Left and Right Dorsolateral Prefrontal Cortex in Major Depression Is Linked to Negative Emotional Judgment: An fMRI Study in Severe Major Depressive Disorder. *Biol. Psychiatry* 63, 369–376.
- Lang, N., Siebner, H. R., Ernst, D., Nitsche, M. A., Paulus, W., Lemon, R. N. & Rothwell, J. C. (2004) Preconditioning with transcranial direct current stimulation sensitizes the motor cortex to rapid-rate transcranial magnetic stimulation and controls the direction of after-effects, *Biological Psychiatry*, 56, 9, 634-639.
- Liu, S., Sheng, J., Li, B., & Zhang, X. (2017). Recent Advances in Non-invasive Brain Stimulation for Major Depressive Disorder. *Frontiers in human neuroscience*, 11, 526.
- Mutz, J. et al. (2018) Efficacy and acceptability of non-invasive brain stimulation for the treatment of adult unipolar and bipolar depression: A systematic review and meta-analysis of randomised sham-controlled trials. *Neuroscience & Biobehavioral Reviews*. 92, 291-303.
- National Institute for Health and Care Excellence (United Kingdom) (2015a). Interventional procedure overview of transcranial direct current stimulation (tDCS) for depression. *NICE Interv. Proced. Program*. 1–42.
- National Institute for Health and Care Excellence (United Kingdom) (2015b). Transcranial direct current stimulation (tDCS) for depression. *NICE Interv. Proced. Guid*. 1–9.
- Palm, U., Schiller, C., Fintescu, Z., Obermeier, M., Keeser, D., Reisinger, E., Pogarell, O., Nitsche, M. A., Möller, H. J., Padberg, F. (2011) Transcranial direct current stimulation in treatment resistant depression: a randomized double-blind, placebo-controlled study. *Brain Stimul*. 5, 3, 242-251.
- Regner, G. G., Pereira, P., Leffa, D. T., de Oliveira, C., Vercelino, R., Fregni, F., & Torres Iraci L. S. (2018) Preclinical to Clinical Translation of Studies of

Transcranial Direct-Current Stimulation in the Treatment of Epilepsy: A Systematic Review, *Frontiers in Neuroscience*, 12, 189.

Roncero, C., Mardigyan, V., Service, E., Singerman, J., Whittaker, K. C., Friedman, M., & Chertkow, H. (2020). Investigation into the effect of transcranial direct current stimulation on cardiac pacemakers, *Brain Stimulation*, 13, 1, 89-95.

Russo, C. et al. (2017). Safety Review of Transcranial Direct Current Stimulation in Stroke. *Neuromodulation: journal of the International Neuromodulation Society*, 20, 3, 215–222.

Sampaio-Junior B, Tortella G, Borrione L, et al. Efficacy and Safety of Transcranial Direct Current Stimulation as an Add-on Treatment for Bipolar Depression: A Randomized Clinical Trial. *JAMA Psychiatry*. 2018;75(2):158–166.

San-Juan, D., Morales-Quezada, L., Orozco Garduño, A. J., Alonso-Vanegas, M., Fernández González-Aragón, M., Espinoza López, D. A., Vázquez Gregorio, R., Anshel, D. J., & Fregni, F. (2015) Transcranial Direct Current Stimulation in Epilepsy, *Brain Stimulation*, 8, 3, 455-464.

Segrave, R. A., Arnold, S., Hoy, K., & Fitzgerald, P. B. (2014). Concurrent cognitive control training augments the antidepressant efficacy of tDCS: a pilot study. *Brain stimulation*, 7, 2, 325–331.

Segrave, R. et al. (2014) Concurrent cognitive control training augments the antidepressant efficacy of tDCS: a pilot study. *Brain Stimul.* 7, 325–33.

Sela, T. & Davidor, M. (2014) Chapter 11: High-Level Cognitive Functions in Healthy Subjects. *The Stimulated Brain, Cognitive Enhancement Using Non-Invasive Brain Stimulation*, 299-329.

Sudbrack-Oliveira, P., Zanichelli Barbosa, M., Thome-Souza, S., Razza, L. B., Gallucci-Neto, J., da Costa Lane Valiengo, L., Brunoni, A. R. (2012). Transcranial direct current stimulation (tDCS) in the management of epilepsy: A systematic review, *Seizure*, 86, 85-95.

Taylor, R., Galvez, V., & Loo, C. (2018). Transcranial magnetic stimulation (TMS) safety: a practical guide for psychiatrists. *Australasian Psychiatry*, 26, 2, 189–192.

Tortella, G. et al. (2015) Transcranial direct current stimulation in psychiatric disorders. *World J. Psychiatry* 5, 88–102.

Valiengo, L., Benseñor, I. M., Goulart, A. C., de Oliveira, J. F., Zanao, T. A., Boggio, P. S., Lotufo, P. A., Fregni, F., Brunoni, A. R. (2013) The sertraline versus electrical current therapy for treating depression clinical study (select-TDCS): results of the crossover and follow-up phases. *Depress Anxiety*. 30, 7, 646-53.

Sooma Oy

Atomitie 5C , FI-00370 Helsinki, Finland

Tel. +358 10 328 9811

Email: info@soomamedical.com

www.soomamedical.com